

## Evategrel

## Chemical Properties


CAS No. : 2760609-74-1

Formula: C<sub>21</sub>H<sub>26</sub>ClNO<sub>7</sub>S

Molecular Weight: 471.95

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Evategrel is an inhibitor of platelet aggregation.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1189 mL	10.5943 mL	21.1887 mL
5 mM	0.4238 mL	2.1189 mL	4.2377 mL
10 mM	0.2119 mL	1.0594 mL	2.1189 mL
50 mM	0.0424 mL	0.2119 mL	0.4238 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481