


Longipenamide B

Chemical Properties

CAS No. : 1383715-05-6
 Formula: C₁₄H₂₅NO₃
 Molecular Weight: 255.35
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Longipenamide B is an alkaloid compound.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.9162 mL	19.581 mL	39.1619 mL
5 mM	0.7832 mL	3.9162 mL	7.8324 mL
10 mM	0.3916 mL	1.9581 mL	3.9162 mL
50 mM	0.0783 mL	0.3916 mL	0.7832 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481