


## Sorocein H

## Chemical Properties

CAS No. : 157382-82-6  
 Formula: C45H44O12  
 Molecular Weight: 776.82  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Sorocein H is a phenolic compound.
-------------	------------------------------------

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2873 mL	6.4365 mL	12.873 mL
5 mM	0.2575 mL	1.2873 mL	2.5746 mL
10 mM	0.1287 mL	0.6436 mL	1.2873 mL
50 mM	0.0257 mL	0.1287 mL	0.2575 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481