


## 3-Methyl-2-butenyl caffeate

## Chemical Properties

CAS No. : 118971-61-2  
 Formula: C<sub>14</sub>H<sub>16</sub>O<sub>4</sub>  
 Molecular Weight: 248.27  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	3-Methyl-2-butenyl caffeate is a type of phenylpropanoid compound.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.0279 mL	20.1394 mL	40.2787 mL
5 mM	0.8056 mL	4.0279 mL	8.0557 mL
10 mM	0.4028 mL	2.0139 mL	4.0279 mL
50 mM	0.0806 mL	0.4028 mL	0.8056 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481