


Schidigerasaponin A3

Chemical Properties

CAS No. :	266997-31-3
Formula:	C45H72O18
Molecular Weight:	901.053
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Schidigerasaponin A3 is a useful organic compound for research related to life sciences. The catalog number is T125158 and the CAS number is 266997-31-3.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1098 mL	5.5491 mL	11.0982 mL
5 mM	0.222 mL	1.1098 mL	2.2196 mL
10 mM	0.111 mL	0.5549 mL	1.1098 mL
50 mM	0.0222 mL	0.111 mL	0.222 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481