

Anti-Sclerostin Monoclonal Antibody

Product Details

Ig Type: Rabbit monoclonal IgG
Reactivity: Human
Conjugation: Unconjugated
Molecular Weight: 150 kDa

Applications

Application: FCM
Recommended: 0.1-0.2 $\mu\text{g}/10\text{E}6$ cells for FCM; 1 ng/ μl for ELISA

Properties

Purity: > 95% as determined by SDS-PAGE.
Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: Sclerostin
Antigen Species: Human
Gene ID: 50964
Biology Area: Drug Target

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481