

## Anti-SLC25A30 Monoclonal Antibody-Biotin

### Product Details

Ig Type: Rabbit monoclonal IgG  
Reactivity: Human  
Conjugation: Biotin  
Molecular Weight: 150 kDa

### Applications

Application: FCM  
Recommended: 0.1-0.2  $\mu\text{g}/10\text{E}6$  cells for FCM; 1 ng/ $\mu\text{l}$  for ELISA

### Properties

Purity: > 95% as determined by SDS-PAGE.  
Stability & Storage: Store at  $-20^{\circ}\text{C}$  or  $-80^{\circ}\text{C}$  for 12 months. Avoid repeated freeze-thaw cycles.  
Shipping: Shipping with blue ice.

### Antigen Details

Immunogen: SLC25A30  
Antigen Species: Human  
Gene ID: 253512  
Biology Area: Other Functions

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481