

Human CD4+ T Cell Isolation Kit (negative selection)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Store at 4°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	TargetMol's Human CD4+ T Cell Isolation Kit (negative selection) provides superparamagnetic beads for the isolation of CD4+ T cells from human peripheral blood mononuclear cells (PBMC) using a negative sorting method. The principle is to use biotin-labelled monoclonal antibody to label non-target cells (non-CD4+ T cells), and then remove the non-target cells by streptavidin-labelled magnetic beads for the purpose of human peripheral blood CD4+ T cell sorting.
-------------	---

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481