

1-Phenyl-1-pyridin-2-ylmethanamine dihydrochloride

Chemical Properties

CAS No. : 39930-11-5


Formula: C₁₂H₁₂N₂

Molecular Weight: 184.242

Store at low temperature

Storage: Pure form: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	1-Phenyl-1-pyridin-2-ylmethanamine dihydrochloride, with CAS No. 39930-11-5, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 1-Phenyl-1-pyridin-2-ylmethanamine dihydrochloride provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.4277 mL	27.1385 mL	54.277 mL
5 mM	1.0855 mL	5.4277 mL	10.8554 mL
10 mM	0.5428 mL	2.7139 mL	5.4277 mL
50 mM	0.1086 mL	0.5428 mL	1.0855 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481