


Catharticin

Chemical Properties

CAS No. :	39723-40-5
Formula:	C ₃₄ H ₄₂ O ₁₉
Molecular Weight:	754.69
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Catharticin is a natural product for research related to life sciences. The catalog number is TN3611 and the CAS number is 39723-40-5.
Targets(IC50)	Others,CDK,NO Synthase,ROS,Tyrosinase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.325 mL	6.625 mL	13.2505 mL
5 mM	0.265 mL	1.325 mL	2.6501 mL
10 mM	0.1325 mL	0.6625 mL	1.325 mL
50 mM	0.0265 mL	0.1325 mL	0.265 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481