

PGLYRP1 Protein, Human, Recombinant (His), Biotinylated

General Information

Synonyms:	UNQ639/PRO1269;TNFSF3L;PGRPS;TAG7;PGRP;PGRP-S;SBBI68;PGLYRP
Protein Construction:	A DNA sequence encoding the Human PGLYRP1 (O75594) (Met1-Pro196) was expressed with a polyhistidine tag at the C-terminus.
Species:	Human
Expression Host:	HEK293 Cells
Accession:	O75594
Molecular Weight:	20.88 kDa (predicted); 24.3 kDa (reducing condition)

QC Testing

Biological Activity:	Activity testing is in progress. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	≥ 95% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	Lyophilized from sterile PBS, pH 7.4. Please contact us for any concerns or special requirements. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween 80 are added as protectants before lyophilization. Please refer to the specific buffer information in the hardcopy of datasheet or the lot-specific COA.

Preparation and Storage

Reconstitution:

Please refer to the lot-specific COA.

Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Peptidoglycan recognition protein 1, also known as Peptidoglycan recognition protein short, PGRP-S, PGLYRP1, PGLYRP, PGRP and TNFSF3L, is a secreted protein that belongs to the N-acetylmuramoyl-L-alanine amidase 2 family. PGLYRP1 / PGLYRP is highly expressed in bone marrow. It is weakly expressed in kidney, liver, small intestine, spleen, thymus, peripheral leukocyte, lung, fetal spleen and neutrophils. PGLYRP1 / PGLYRP is a pattern

receptor that binds to murein peptidoglycans (PGN) of Gram-positive bacteria. It has bactericidal activity towards Gram-positive bacteria. PGLYRP1 / PGLYRP may kill Gram-positive bacteria by interfering with peptidoglycan biosynthesis. It binds also to Gram-negative bacteria, and has bacteriostatic activity towards Gram-negative bacteria. Peptidoglycan recognition proteins (PGRPs or PGLYRPs) are innate immunity proteins that are conserved from insects to mammals, recognize bacterial peptidoglycan, and function in antibacterial immunity and inflammation. Mammals have four PGRPs: PGLYRP1, PGLYRP2, PGLYRP3, and PGLYRP4. They are secreted proteins expressed in polymorphonuclear leukocytes (PGLYRP1), liver (PGLYRP2), or on body surfaces, mucous membranes, and in secretions (saliva, sweat) (PGLYRP3 and PGLYRP4). All PGRPs recognize bacterial peptidoglycan. The PGRPs likely play a role both in antibacterial defenses and several inflammatory diseases. They modulate local inflammatory responses in tissues (such as arthritic joints) and there is evidence for association of PGRPs with inflammatory diseases, such as psoriasis.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481