

Apolipoprotein E/APOE3 Protein, Human, Recombinant (Isoform E4, C130R)

General Information

Synonyms:	LDLCQ5;AD2;LPG;APOE;APO-E
Protein Construction:	A DNA sequence encoding the Human ApoE (P02649) (Met1-His317, C130R) was expressed.
Species:	Human
Expression Host:	HEK293 Cells
Accession:	P02649
Molecular Weight:	34.3 kDa (predicted); 34 kDa (reducing condition)

QC Testing

Biological Activity:	Activity testing is in progress. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	≥ 90% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	Lyophilized from sterile PBS, pH 7.4. Please contact us for any concerns or special requirements. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween 80 are added as protectants before lyophilization. Please refer to the specific buffer information in the hardcopy of datasheet or the lot-specific COA.

Preparation and Storage

Reconstitution:

Please refer to the lot-specific COA.

Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Apolipoprotein E (ApoE) is a 34.2 kDa glycosylated protein with 299 amino acid residues. There are three isoforms in human (apoE2, apoE3, and apoE4) due to different amino acid residues at positions 112 and 158. ApoE is synthesized predominantly in the liver, but also by cells in the spleen, brain, lung, kidney, ovary, adrenal, and muscle tissues. Hepatic parenchyma cells are the main apoE producing cells in mammalian body, probably accounting for two thirds to three fourths of the plasma apoE. In the nervous system, apoE mRNA is present in

neurons, astrocytes, ependymal cells, nonmyelinating Schwann cells, but not in microglia, oligodendroglia, choroidal cells, or myelinating Schwann cells. ApoE produced by mammalian cells exists in different forms, monomers, dimers, modified, unmodified, lipid-rich, and lipid-poor, and so forth. ApoE plays a double-role in immune responses. Both apoE containing lipoproteins and multimers of synthetic apoE peptides inhibited proliferation of cultured lymphocytes by inhibiting DNA synthesis and reducing phospholipid turnover in T cells. ApoE can also affect innate and acquired immune responses in vitro by its ability to suppress stimulation of cultured neutrophils. ApoE can bind lipopolysaccharide (LPS), attenuate the inflammatory response, and thus reduce LPS induced lethality. Injection of LPS stimulated higher expression of inflammatory cytokines like interleukin (IL)-1 β , IL-12, and interferon- γ (IFN- γ), as well as IL-6.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481