

GIPR Protein, Human, Recombinant (HEK293, His)

General Information

Synonyms:	PGQTL2
Protein Construction:	A DNA sequence encoding the Human GIPR (P48546-1) (Gly26-Gln138) was expressed with a polyhistidine tag at the C-terminus.
Species:	Human
Expression Host:	HEK293 Cells
Accession:	P48546-1
Molecular Weight:	14.42 kDa (predicted); 25.7 kDa (reducing condition, due to glycosylation)

QC Testing

Biological Activity:	Immobilized anti-GIPR, Human IgG1 at 2 µg/mL (100 µL/well) can bind GIPR Protein, Human, Recombinant (HEK293, His) (Cat#TMPY-07153), the EC50 is 1.5-4.5 ng/mL.
Purity:	≥ 95% as determined by SDS-PAGE. ≥ 95% as determined by SEC-HPLC.
Endotoxin:	< 1.0 EU/µg of the protein as determined by the LAL method.
Formulation:	Lyophilized from sterile PBS, pH 7.4. Please contact us for any concerns or special requirements. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween 80 are added as protectants before lyophilization. Please refer to the specific buffer information in the hardcopy of datasheet or the lot-specific COA.

Preparation and Storage

Reconstitution:

Please refer to the lot-specific COA.

Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481