

Compound T130436

Chemical Properties


CAS No. :

Formula: C₂₉H₃₄O₁₇

Molecular Weight: 654.57

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	6-hydroxy-3-(3-hydroxy-4-([3,4,5-trihydroxy-6-((3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy)methyl)oxan-2-yl]oxy)phenyl)-5,7-dimethoxy-4H-chromen-4-one is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5277 mL	7.6386 mL	15.2772 mL
5 mM	0.3055 mL	1.5277 mL	3.0554 mL
10 mM	0.1528 mL	0.7639 mL	1.5277 mL
50 mM	0.0306 mL	0.1528 mL	0.3055 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481