


## Photoregulin3

## Chemical Properties

CAS No. :	785708-33-0
Formula:	C <sub>22</sub> H <sub>23</sub> N <sub>3</sub> O <sub>3</sub>
Molecular Weight:	377.44
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Photoregulin3 is a inhibitor of rod photoreceptor gene expression, potentially through Nr2e3 modulation.
Targets(IC50)	Others,DNA/RNA Synthesis

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6494 mL	13.2471 mL	26.4943 mL
5 mM	0.5299 mL	2.6494 mL	5.2989 mL
10 mM	0.2649 mL	1.3247 mL	2.6494 mL
50 mM	0.053 mL	0.2649 mL	0.5299 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481