

E3 Ligase Ligand-linker Conjugate 136

Chemical Properties

CAS No. : 2761394-06-1


Formula: C₂₂H₂₂N₄O₄

Molecular Weight: 406.44

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	E3 Ligase Ligand-linker Conjugate 136 is an E3 ubiquitin ligase ligand-linker conjugate (E3 Ligase Ligand-linker conjugate) used for the synthesis of c-Myc inhibitor7.
-------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4604 mL	12.3019 mL	24.6039 mL
5 mM	0.4921 mL	2.4604 mL	4.9208 mL
10 mM	0.246 mL	1.2302 mL	2.4604 mL
50 mM	0.0492 mL	0.246 mL	0.4921 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481