

## TRIB2 ligand-1

## Chemical Properties

CAS No. : 2652012-89-8


Formula: C<sub>22</sub>H<sub>18</sub>FN<sub>5</sub>O

Molecular Weight: 387.41

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	TRIB2 ligand-1 is a ligand for TRIB2 (Ligands for Target Protein for PROTACs). This compound serves as a precursor for the synthesis of PROTAC TRIB2 degrader-1.
Targets(IC50)	Ligands for Target Protein for PROTAC

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5812 mL	12.9062 mL	25.8124 mL
5 mM	0.5162 mL	2.5812 mL	5.1625 mL
10 mM	0.2581 mL	1.2906 mL	2.5812 mL
50 mM	0.0516 mL	0.2581 mL	0.5162 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481