

## Vislarafusp alfa

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

-20°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Vislarafusp alfa is a humanized IgG1k antibody that targets EGFR/CD47, with HumanIgG1kappa, Isotype Control serving as its corresponding isotype control.
Targets(IC50)	EGFR

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481