

## p62-ZZ Ligand-Linker Conjugate 1

### Chemical Properties

CAS No. : 2410081-73-9


Formula: C<sub>31</sub>H<sub>42</sub>N<sub>2</sub>O<sub>6</sub>

Molecular Weight: 538.68

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



### Biological Description

Description	p62-ZZ Ligand-Linker Conjugate 1 is a conjugate of the p62-ZZ domain ligand and linker. This compound is utilized in the synthesis of AUTOTACVinclozolinM2-2204.
Targets(IC50)	Autophagy

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8564 mL	9.2819 mL	18.5639 mL
5 mM	0.3713 mL	1.8564 mL	3.7128 mL
10 mM	0.1856 mL	0.9282 mL	1.8564 mL
50 mM	0.0371 mL	0.1856 mL	0.3713 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481