

E3 Ligase Ligand-linker Conjugate 137

Chemical Properties

CAS No. :

Formula: C₂₅H₃₃N₅O₆

Molecular Weight: 499.56

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	E3 Ligase Ligand-linker Conjugate 137 serves as both a linker and an E3 ligase ligand for PROTAC METTL3-14 degrader 1, employed in the study of acute myeloid leukemia.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0018 mL	10.0088 mL	20.0176 mL
5 mM	0.4004 mL	2.0018 mL	4.0035 mL
10 mM	0.2002 mL	1.0009 mL	2.0018 mL
50 mM	0.040 mL	0.2002 mL	0.4004 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481