

MTH1 ligand 1

Chemical Properties

CAS No. : 2412986-35-5


Formula: C₂₃H₁₈N₄O₃

Molecular Weight: 398.41

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	MTH1 ligand 1 functions as a target protein ligand for MTH1, and is utilized in the synthesis of PROTACaTAG 2139.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates,DNA/RNA Synthesis,Ligands for Target Protein for PROTAC

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.510 mL	12.5499 mL	25.0998 mL
5 mM	0.502 mL	2.510 mL	5.020 mL
10 mM	0.251 mL	1.255 mL	2.510 mL
50 mM	0.0502 mL	0.251 mL	0.502 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481