

## Avitotamig

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Avitotamig is a bispecific antibody of the [H- $\gamma$ 1 _ L- $\kappa$ -scFv heavy - $\kappa$ ]-dimer type that targets CD3E/CD33.
Targets(IC50)	Integrin

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481