


PROTAC BRM/BRG1 degrader-2

Chemical Properties

CAS No. :	2933124-32-2
Formula:	C50H60N12O5S
Molecular Weight:	941.15
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	PROTAC BRM/BRG1 degrader-2 (Example.004) is a potent PROTAC BRM/BRG1 degrading agent exhibiting a DC50 of less than 10 nM and an IC50 of 15 nM in A549 cells.
Targets(IC50)	Epigenetic Reader Domain,PROTACs

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0625 mL	5.3126 mL	10.6253 mL
5 mM	0.2125 mL	1.0625 mL	2.1251 mL
10 mM	0.1063 mL	0.5313 mL	1.0625 mL
50 mM	0.0213 mL	0.1063 mL	0.2125 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481