


E3 ligase Ligand 43

Chemical Properties

CAS No. :	2821795-71-3
Formula:	C ₂₂ H ₂₈ N ₆ O ₄ S
Molecular Weight:	472.56
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	E3 Ligase Ligand 43 is an E3 ubiquitin ligase ligand for ACB13, used in the synthesis of PROTACs.
-------------	---------------------------------------------------------------------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1161 mL	10.5807 mL	21.1613 mL
5 mM	0.4232 mL	2.1161 mL	4.2323 mL
10 mM	0.2116 mL	1.0581 mL	2.1161 mL
50 mM	0.0423 mL	0.2116 mL	0.4232 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481