

## SMARCA2 ligand-8

## Chemical Properties

CAS No. : 2933254-33-0


Formula: C<sub>12</sub>H<sub>9</sub>N<sub>4</sub>O

Molecular Weight: 352.13

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	SMARCA2 ligand-8 serves as a ligand for the target protein SMARCA2, utilized in the synthesis of PROTAC SMARCA2/4-degrader-35.
Targets(IC50)	Ligands for Target Protein for PROTAC

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8399 mL	14.1993 mL	28.3986 mL
5 mM	0.568 mL	2.8399 mL	5.6797 mL
10 mM	0.284 mL	1.4199 mL	2.8399 mL
50 mM	0.0568 mL	0.284 mL	0.568 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481