

Obrixtamig

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Obrixtamig is a chimeric human antibody of the L-k-G1-h-CH2-CH3_L-λ-G1-h-CH2-CH3 type, targeting CD3E/HEK.
Targets(IC50)	Integrin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481