

## PROTAC SMARCA2 degrader-29

## Chemical Properties

CAS No. : 2865193-25-3


Formula: C47H54N8O5S2

Molecular Weight: 875.11

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	PROTAC SMARCA2 degrader-29 (Example 87) is a potent PROTAC SMARCA2 (BRM) degrading agent with a DC50 ranging between 10-100 nM (BRM) and over 1 μM (BRG1, also known as SMARCA2).
Targets(IC50)	Epigenetic Reader Domain, Neurotensin Receptor, PROTACs

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1427 mL	5.7136 mL	11.4271 mL
5 mM	0.2285 mL	1.1427 mL	2.2854 mL
10 mM	0.1143 mL	0.5714 mL	1.1427 mL
50 mM	0.0229 mL	0.1143 mL	0.2285 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481