


PI3Ky ligand 1

Chemical Properties

CAS No. :	2769753-45-7
Formula:	C ₂₆ H ₂₉ N ₅ O ₃ S
Molecular Weight:	491.61
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	PI3Ky ligand 1 serves as a PROTAC target protein ligand (Ligand for Target Protein for PROTAC).
Targets(IC50)	Ligands for Target Protein for PROTAC,PI3K

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0341 mL	10.1707 mL	20.3413 mL
5 mM	0.4068 mL	2.0341 mL	4.0683 mL
10 mM	0.2034 mL	1.0171 mL	2.0341 mL
50 mM	0.0407 mL	0.2034 mL	0.4068 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481