


## PROTAC IRAK4 ligand-5

## Chemical Properties

CAS No. :	2654056-24-1
Formula:	C <sub>24</sub> H <sub>24</sub> F <sub>3</sub> N <sub>3</sub> O <sub>3</sub> S
Molecular Weight:	491.53
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	PROTAC IRAK4 ligand-5, a ligand for the target protein of PROTAC, is instrumental in the synthesis of KT-413.
Targets(IC50)	Ligands for Target Protein for PROTAC

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0345 mL	10.1723 mL	20.3446 mL
5 mM	0.4069 mL	2.0345 mL	4.0689 mL
10 mM	0.2034 mL	1.0172 mL	2.0345 mL
50 mM	0.0407 mL	0.2034 mL	0.4069 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481