

Me-Thalidomide-O-C8-NH2

Chemical Properties

CAS No. : 2550393-20-7


Formula: C22H29N3O5

Molecular Weight: 415.48

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Me-Thalidomide-O-C8-NH2 is an E3 ubiquitin ligase ligand-linker conjugate (E3 Ligase Ligand-linker conjugate). This compound is utilized in the synthesis of MS154N.
Targets(IC50)	Antibacterial, E3 Ligase Ligand-Linker Conjugates

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4069 mL	12.0343 mL	24.0685 mL
5 mM	0.4814 mL	2.4069 mL	4.8137 mL
10 mM	0.2407 mL	1.2034 mL	2.4069 mL
50 mM	0.0481 mL	0.2407 mL	0.4814 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481