

Emunkitug

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Emunkitug is a humanized IgG1κ monoclonal antibody that targets TNFRSF1B. Its corresponding isotype control is HumanIgG1kappa, Isotype Control.
Targets(IC50)	TNF

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481