


BRM/BRG1 ligand 1

Chemical Properties

CAS No. :	1893977-70-2
Formula:	C ₁₀ H ₈ BrN ₃ O
Molecular Weight:	266.09
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	BRM/BRG1 Ligand 1 is a BRM/BRG1 ligand utilized for the synthesis of PROTAC BRM/BRG1 degrader-1.
Targets(IC50)	Epigenetic Reader Domain,Ligands for Target Protein for PROTAC

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.7581 mL	18.7906 mL	37.5813 mL
5 mM	0.7516 mL	3.7581 mL	7.5163 mL
10 mM	0.3758 mL	1.8791 mL	3.7581 mL
50 mM	0.0752 mL	0.3758 mL	0.7516 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481