


AR ligand-30

Chemical Properties

CAS No. :	2505498-73-5
Formula:	C ₁₈ H ₁₆ Cl ₂ N ₄ O ₂
Molecular Weight:	391.25
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	AR ligand-30 serves as the target protein ligand for PROTAC Bavdegalutamide, utilized in prostate cancer research.
Targets(IC50)	Antifungal, Ligands for Target Protein for PROTAC

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5559 mL	12.7796 mL	25.5591 mL
5 mM	0.5112 mL	2.5559 mL	5.1118 mL
10 mM	0.2556 mL	1.278 mL	2.5559 mL
50 mM	0.0511 mL	0.2556 mL	0.5112 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481