

Lenalidomide-ethyne-C2-Pip-CO-Pip

Chemical Properties

CAS No. :

Formula: C₂₆H₃₂N₆O₄

Molecular Weight: 492.57

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Lenalidomide-ethyne-C2-Pip-CO-Pip (intermediate 3) is an E3 ubiquitin ligase ligand-linker conjugate. It serves as a precursor in the synthesis of QA-68.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates,Molecular Glues

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0302 mL	10.1508 mL	20.3017 mL
5 mM	0.406 mL	2.0302 mL	4.0603 mL
10 mM	0.203 mL	1.0151 mL	2.0302 mL
50 mM	0.0406 mL	0.203 mL	0.406 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481