

ASP-8731

Chemical Properties


CAS No. : 2488255-42-9

Formula: C₂₀H₂₁N₅O₄

Molecular Weight: 395.41

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	ASP8731 is a selective inhibitor of BACH1 that prevents inflammation and vascular occlusion, and induces fetal hemoglobin in sickle cell disease.
Targets(IC50)	Others,ROS

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.529 mL	12.6451 mL	25.2902 mL
5 mM	0.5058 mL	2.529 mL	5.058 mL
10 mM	0.2529 mL	1.2645 mL	2.529 mL
50 mM	0.0506 mL	0.2529 mL	0.5058 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481