

E3 Ligase Ligand-linker Conjugate 125

Chemical Properties

CAS No. : 2918284-14-5


Formula: C₂₂H₂₃N₅O₄

Molecular Weight: 421.45

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

E3 Ligase Ligand-linker Conjugate 125 is an E3 ubiquitin ligase ligand-linker conjugate (E3 Ligase Ligand-linker Conjugate) used for synthesizing PROTAC IRAK4degrader-12.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3728 mL	11.8638 mL	23.7276 mL
5 mM	0.4746 mL	2.3728 mL	4.7455 mL
10 mM	0.2373 mL	1.1864 mL	2.3728 mL
50 mM	0.0475 mL	0.2373 mL	0.4746 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481