

## PROTAC SMARCA2 degrader-28

## Chemical Properties


CAS No. : 2641510-88-3

Formula: C<sub>50</sub>H<sub>51</sub>F<sub>7</sub>N<sub>10</sub>O<sub>10</sub>

Molecular Weight: 1084.99

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	PROTAC SMARCA2 degrader-28 (Compound 158) is a potent SMARCA2 PROTAC degrader that achieves degradation of SMARCA2 in HiBiT A549 cells, exhibiting a DC50 of 3 nM.
Targets(IC50)	Epigenetic Reader Domain,PROTACs

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9217 mL	4.6083 mL	9.2167 mL
5 mM	0.1843 mL	0.9217 mL	1.8433 mL
10 mM	0.0922 mL	0.4608 mL	0.9217 mL
50 mM	0.0184 mL	0.0922 mL	0.1843 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481