

## Amtabafusp alfa

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

|               |  |
|---------------|--|
| Description   | Amtabafusp alfa is a humanized IgG1λ1 antibody targeting CD3E/CD4. Its isotype control corresponds to HumanIgG1lambda1, Isotype Control. |
| Targets(IC50) | Integrin   |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481