

ADC Control Human IgG1-DM4

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	ADC Control Human IgG1-DM4 serves as the isotype control for ADC Human IgG1-DM4. It comprises an antibody portion identified as HumanIgG1kappa, Isotype Control, whereas the toxic molecule and linker components are SPDB-DM4.
Targets(IC50)	Antibody-Drug Conjugates (ADCs)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481