

PROTAC BTK Degradator-11

Chemical Properties

CAS No. : 2736508-94-2


Formula: C₄₈H₅₅N₁₁O₄

Molecular Weight: 850.02

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	PROTAC BTK Degradator-11 is a PROTAC degrader that targets and degrades BTK with a DC50 of 1.7 nM. It is utilized in oncological research.
Targets(IC50)	BTK,PROTACs

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1764 mL	5.8822 mL	11.7644 mL
5 mM	0.2353 mL	1.1764 mL	2.3529 mL
10 mM	0.1176 mL	0.5882 mL	1.1764 mL
50 mM	0.0235 mL	0.1176 mL	0.2353 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481