

## (3S)Lenalidomide-piperazine-C-piperidine

## Chemical Properties

CAS No. : 2520107-24-6


Formula: C23H31N5O3

Molecular Weight: 425.52

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(3S)Lenalidomide-piperazine-C-piperidine is an E3 ubiquitin ligase ligand-linker conjugate (E3 Ligase Ligand-Linker Conjugates). This compound can be utilized in the synthesis of PROTAC ER Degradar-11.
Targets(IC50)	Antibacterial, Antibiotic

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3501 mL	11.7503 mL	23.5007 mL
5 mM	0.470 mL	2.3501 mL	4.7001 mL
10 mM	0.235 mL	1.175 mL	2.3501 mL
50 mM	0.047 mL	0.235 mL	0.470 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481