

BRD9 ligand-6

Chemical Properties

CAS No. : 2460248-49-9


Formula: C₂₁H₂₁N₃O₄

Molecular Weight: 379.41

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	BRD9 ligand-6 serves as a ligand for the target protein in PROTAC applications. It is utilized in the synthesis of FHD-609.
Targets(IC50)	5-HT Receptor, Ligands for Target Protein for PROTAC

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6357 mL	13.1784 mL	26.3567 mL
5 mM	0.5271 mL	2.6357 mL	5.2713 mL
10 mM	0.2636 mL	1.3178 mL	2.6357 mL
50 mM	0.0527 mL	0.2636 mL	0.5271 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481