


## E3 Ligase Ligand-linker Conjugate 140

### Chemical Properties

CAS No. :	2753651-22-6
Formula:	C30H34N6O5
Molecular Weight:	558.63
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



### Biological Description

Description	E3 Ligase Ligand-linker Conjugate 140 (intermediate I-11) is an E3 ubiquitin ligase ligand-linker conjugate used for synthesizing PROTACAR Degrader-6.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7901 mL	8.9505 mL	17.9009 mL
5 mM	0.358 mL	1.7901 mL	3.5802 mL
10 mM	0.179 mL	0.895 mL	1.7901 mL
50 mM	0.0358 mL	0.179 mL	0.358 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481