


## Tipindole

## Chemical Properties

CAS No. : 7489-66-9  
 Formula: C<sub>16</sub>H<sub>20</sub>N<sub>2</sub>O<sub>2</sub>S  
 Molecular Weight: 304.41  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Tipindole is a serotonin antagonist utilized in research related to depression.
Targets(IC50)	5-HT Receptor

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.285 mL	16.4252 mL	32.8504 mL
5 mM	0.657 mL	3.285 mL	6.5701 mL
10 mM	0.3285 mL	1.6425 mL	3.285 mL
50 mM	0.0657 mL	0.3285 mL	0.657 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481