

CDK1-IN-6

## Chemical Properties


CAS No. :

Formula: C<sub>21</sub>H<sub>22</sub>N<sub>4</sub>O

Molecular Weight: 346.43

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	CDK1-IN-6 (Ligand 3) is an effective inhibitor of CDK1 and shows potential for use in cancer research.
Targets(IC50)	CDK

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8866 mL	14.4329 mL	28.8659 mL
5 mM	0.5773 mL	2.8866 mL	5.7732 mL
10 mM	0.2887 mL	1.4433 mL	2.8866 mL
50 mM	0.0577 mL	0.2887 mL	0.5773 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481