

Sulfo SMCC R-phycoerythrin

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	SulfoSMCCR-phycoerythrin is a conjugate composed of the protein crosslinker SMCC and R-PE (R-Phycoerythrin), designed for labeling proteins to carry red fluorescence. SMCC functions to bind antigen-coupled spleen cells, thereby inducing antigen-specific immune responses.
Targets(IC50)	Epigenetic Reader Domain

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481