

Veбанvibart

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Veбанvibart is a humanized IgG1k monoclonal antibody targeting SARS-CoV-2, with HumanIgG1kappa, Isotype Control serving as its corresponding isotype control.
Targets(IC50)	SARS-CoV

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481