

Zaltenibart

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Zaltenibart is a chimeric humanized antibody of the IgG4k type, specifically targeting MASP3. The corresponding isotype control for this antibody is: HumanIgG4(S228P) kappa, Isotype Control.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481