

BNN-20

## Chemical Properties


CAS No. : 847-75-6

Formula: C<sub>20</sub>H<sub>30</sub>O<sub>2</sub>

Molecular Weight: 302.45

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	BNN-20 is a synthetic micro-neurotrophic factor that mimics Brain-Derived Neurotrophic Factor (BDNF) and exhibits region-specific neuroprotective and neurogenesis-promoting effects. It is applicable in research related to Parkinson's disease.
Targets(IC50)	Trk receptor

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3063 mL	16.5317 mL	33.0633 mL
5 mM	0.6613 mL	3.3063 mL	6.6127 mL
10 mM	0.3306 mL	1.6532 mL	3.3063 mL
50 mM	0.0661 mL	0.3306 mL	0.6613 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481