


## Sarafotoxin S6b

## Chemical Properties

CAS No. :	116303-65-2
Formula:	C110H163N27O34S5
Molecular Weight:	2567.96
Storage:	Keep away from moisture Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	non-selective endothelin receptor agonist
Targets(IC50)	Endothelin Receptor,Others

## Solubility Information

Solubility	20% acetonitrile: 1 mg/mL (0.39 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.3894 mL	1.9471 mL	3.8941 mL
5 mM	0.0779 mL	0.3894 mL	0.7788 mL
10 mM	0.0389 mL	0.1947 mL	0.3894 mL
50 mM	0.0078 mL	0.0389 mL	0.0779 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481